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Quartz Retarders: low order

The retarders are single plates with a path difference between the principal polarization axes equal to the desired wavelength fraction added with an integer multiple of the specified design wavelength λ . The integer part k (order number) is chosen to yield a mechanical thickness of about 0.1 mm to 0.2 mm. The resulting order numbers are $k \geq 5$ at $\lambda=250$ nm and $k \leq 2$ at $\lambda=1000$ nm.

The wavefront distortion is $\lambda/10$ (at $\lambda = 550$ nm) and the wedge error is smaller than $2''$. The accuracy of the path difference is ± 3 nm.

The spectral transmission range of quartz is 180 nm - 2.7 μ m. The usable spectral range for a single low-order plate is much smaller and depends on the required accuracy of the retardation as well as on the possible order number k . At $\lambda=633$ nm ($k=2$) wavelength changes of about ± 1 nm will cause a deviation of the path difference smaller than the production tolerance of ± 3 nm. At $\lambda=355$ nm ($k=4$) a wavelength change of only 0.5 nm will result in the same deviation. Although a deviation of the path difference which is much larger than the production tolerance is acceptable in many applications we recommend the use of zero order plates for wavelengths below 350 nm. Nevertheless, we will deliver low-order plates for short wavelengths on request. Below 230 nm these plates will be manufactured from synthetic quartz crystal.


Owing to the small thickness of our low-order plates, the path difference of the plate is neither very dependent on temperature nor on the aperture angle of the incident beam of light. The deviation ΔR of the path difference from the nominal value R (both in nm) due to a temperature change by ΔT (in K) or due to tilt by $\Delta \varphi$ (in degrees) can be estimated as follows:

$$\Delta R \approx -(0.1 \text{ nm/K}) \cdot \Delta T$$

$$\Delta R \approx \pm(0.1 \text{ nm/deg}^2) \cdot \Delta \varphi^2$$

The deviations will be smaller than the production tolerance for a temperature change up to ± 30 K or for a tilt up to $\pm 5.5^\circ$.

The plates are AR coated with double layers at the design wavelength ($R < 0.1\%$ per surface) and are delivered with a cylindrical mount made of anodized aluminum and engraved with the direction of the fast axis, the path difference and the design wavelength.

Öffnung/Aperture Durchmesser/diameter	Fassung/Holder		$\lambda/2$	$\lambda/4$
	Durchmesser/diameter	Länge/length		
9.5 mm	25 mm	10 mm	RLQ 2.10	RLQ 4.10
14.5 mm	25 mm	10 mm	RLQ 2.15	RLQ 4.15
19.5 mm	30 mm	15 mm	RLQ 2.20	RLQ 4.20
24.5 mm	40 mm	15 mm	RLQ 2.25	RLQ 4.25

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Important note:

Please specify required wavelength when ordering.
The prices apply to standard wavelengths according to the following table.

Standardwellenlängen für low order Platten aus Quarz Standard wavelengths for low order quartz retarders						
355 nm	375 nm	390 nm	395 nm	400 nm	405 nm	425 nm
440 nm	488 nm	514 nm	532 nm	589 nm	633 nm	650 nm
670 nm	685 nm	780 nm	785 nm	790 nm	795 nm	800 nm
810 nm	830 nm	852 nm	946 nm	1030 nm	1047 nm	1053 nm
1064 nm	1083 nm	1300 nm	1320 nm	1540 nm	1550 nm	2023 nm

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Available on special order:

non-standard wavelengths have a surcharge:

clear aperture 9.5 mm and 14.5 mm RLQ 0.1
clear aperture 19.5 mm and 24.5 mm RLQ 0.2

$\lambda/8$ and other retardation values on request
other dimensions on request
retarders for 2 wavelengths on request
accuracy of retardation ± 1 nm RLQ 0.5